

MP-225 MOTORIZED MICROMANIPULATOR

MP-225 电动显微操作器



The **MP-225** is one of the newer products in the Sutter Instrument micromanipulator line. Production and design changes have allowed us to produce a motorized manipulator that is a more affordable alternative to the industry standard **MP-285**. While the **MP-225** feature set is less comprehensive than the **MP-285**, it includes the most popular features with a simpler user interface. The mechanical design utilizes a miniature stepper motor and integral anti-backlash gear head. Pre-loaded ball bearing slides provide smooth movement throughout the 25mm of travel. The controller uses low-noise, linear-drive output circuitry identical to that found in the **MP-285**. The methodology for mounting pipette holders and headstages used with the **MP-285** has been maintained in the **MP-225** to allow for cross compatibility.

The **MP-225** is designed primarily for positioning patch and intracellular recording pipettes. We have retained and refined the features most desired for this type of work. An extended version of the popular rotary optical encoder (ROE) is the sole input device available with the **MP-225**. Like the **MP-285**, the manipulator has a synthetic 4th axis for diagonal advancement of the pipette; 16 different angles are selectable via DIP switches. Speed and resolution of movement are easily selected with a multiple position thumbwheel, allowing fast/coarse movement and slow/ultra-fine movement in 10 increments. Two commonly used robotic movements have been incorporated for user convenience. A single button press can initiate a move to a Home position for pipette exchange or to a user defined Work position for quick location of the pipette

马普科学仪器有限公司

Tel: 020-8767 9617, 8767 9631; Fax: 020-8767 9635
<http://www.mapusci.com> E-mail: info@mapusci.com

near the recording location. A display on the ROE gives position location. As all controls are located on the ROE, the controller can be moved to a less accessible area of your setup and does not need to occupy prime space in an equipment rack.

As always, our technical support team is available to address your particular needs and answer all questions before and after your purchase.

马普科学仪器有限公司

Tel: 020-8767 9617, 8767 9631; Fax: 020-8767 9635
<http://www.mapusci.com> E-mail: info@mapusci.com

FEATURES**> MP-225**

Highly stable for experiments intolerant of pipette drift

Submicron (62.5nm) minimal resolution for fine movement

Convenient thumbwheel selects resolution/speed of movement

25mm of motorized travel on all three axes

4th axis for coaxial movement of pipette, angle selected by DIP switches on ROE

ROE button press actuates move to Home position for pipette exchange

ROE button press actuates move to Work position near recording location

Continuous display (in microns) of axes positions located on ROE

DIP switches on ROE select direction of movement produced by turn of ROE knob

Modularized, compact design easily adaptable to your setup

Universal mounting system for headstage or pipette holder

Mounting adapters included with manipulator

SPECIFICATIONS**> MP-225****Travel**

25mm on all three axes

Resolution

Six microstep sizes selectable ($\mu\text{m}/\text{ustep}$): 0.0625, 0.125, 0.25, 0.5, 1.0 and 2.0

Finer movement settings use the 62.5nm microstep size but fewer microsteps are commanded per encoder knob turn

Maximum Speed

2.0mm/sec

Long Term Stability

1-2 $\mu\text{m}/\text{hour}$ maximum

Drive Mechanism

Integral miniature stepper motor anti-backlash gearhead

Dimensions

Mechanical: 4in x 5.5in x 6in
10cm x 15cm x 15.5cm

Controller: 4in x 16in x 12.25in
10cm x 40.5cm x 31cm

Weight

Manipulator: 2.95lbs /1.3kg

Controller: 10lb 11oz/4.5kg

Electrical

115/230 Volts

50/60 Hertz power line



马普科学仪器有限公司

Tel: 020-8767 9617, 8767 9631; Fax: 020-8767 9635

<http://www.mapusci.com> E-mail: info@mapusci.com

US PRICES**> MP-225**

*International prices are 5% higher.
Prices subject to change without notice.*

Motorized Micromanipulator

MP-225 Includes manipulator, table top controller unit, extended ROE, hinged headstage mount, cables, rod holder, 4in dovetail extension, mounting adapter plate, screws, hex wrench, power cord and manual*

For detailed information on mounting our micromanipulators, refer to the Translators, Stages, & Mounting section or phone for assistance.

** Indicate right- or left-handed set-up when ordering.*

ACCESSORIES

X285204	4 inch dovetail extension
X285210*	Mounting adapter plate
225RBI	Rotating base for MP-225
X285300	Right angle adapter
X285305	Z-axis vertical extension
X285310	Z-axis horizontal extension

** For use with MT and MD stands and platforms, or any surface with one inch centered holes.*

马普科学仪器有限公司

Tel: 020-8767 9617, 8767 9631; Fax: 020-8767 9635
http://www.mapusci.com E-mail: info@mapusci.com